



Instructions for use

Ambu Junior

ENGLISH

1. Introduction

The Ambu Junior is an instruction and training manikin for the simulation of realistic and correct anatomical conditions during cardiopulmonary resuscitation by means of insufflation and external cardiac compression.

The Ambu Junior features a unique, patented hygienic system under which each trainee gets his own face piece and head bag, thereby minimising any risk of cross-infection between trainees.

With the special hygienic system, disassembly and internal cleaning and disinfection are rendered unnecessary both during and after training.

The Ambu Junior provides a lifelike representation of the human anatomy, particularly of those features important to training in modern resuscitation techniques.

2. Hygienic System ^①

The head bag has been mounted in the skull and the mouth/nose piece applied.

Insufflation has just started, see 1.1.

A: Trainee's air **B:** Ambient air

When air is blown through the mouth of the mouth/nose piece, the head bag fills, thereby displacing the air in the skull through the valve, which opens if the head is hyperextended.

Insufflation is concluded and expiration has just started, see 1.2.

NOTE: If a leakage occurs around the mouth/nose piece or if insufflation proves difficult, it is normally due to incorrect positioning when lifting the chin or to insufficient hyperextension of the head.

To comply with the requirements of the hygienic system, and to achieve a proper seal and accurate training results, a head bag must always be used.

All trainees must use their own mouth/nose piece and head bag only during training.

3. Preparation for training.

Fitting of the mouth/nose piece and the head bag ②

Fold the head bag carefully (2.1) and place it in the mouth opening (2.2). Fold down the edge of the head bag and place the sealing ring over it (2.3). Place the mouth/nose piece in the opening of the head bag and press it against the head to obtain a tight seal (2.4).

At the end of the instruction period, discard the head bag and collect the mouth/nose pieces for subsequent cleaning (see chapter 4).

NOTE: Do not twist the bag while folding, as this may cause a considerable increase in resistance to insufflation and expiration.

Removal of mouth/nose piece and head bag ③

Pull the mouth/nose piece upwards. Remove the sealing ring gently (3.1). Remove the head bag slowly upward (3.2). Do not close the head bag with the hand (3.3), as this will make it almost impossible to pull out the head bag.

4. Cleaning and disinfecting

The following procedure should be used:

1. Clean in mild soapy water.
2. Rinse thoroughly in clear water.
3. Disinfect as per C.D.C. protocols.

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Ambu A/S
Simulaids, Inc.